IMPORTANT INFORMATION – PLEASE READ BEFORE USE  
INFORMATIONS IMPORTANTES – À LIRE AVANT UTILISATION  
INFORMAZIONI IMPORTANTI – LEGGERE PRIMA DELL'USO  
WICHTIGE INFORMATIONEN – VOR GEBRAUCH BITTE LESEN  
INFORMACIÓN IMPORTANTE – LÉALA ANTES DE UTILIZAR ESTE PRODUCTO  
重要な情報 – 使用前にお読みください  
ВАЖНА ИНФОРМАЦИЯ – МОЛЯ, ПРОЧЕТЕТЕ ПРЕДИ УПОТРЕБА  
VAŽNE INFORMACIJE – PROČITATI PRIJE UPORABE  
DŮLEŽITÉ INFORMACE – PŘEČTĚTE SI PŘED POUŽITÍM  
VIGTIGE OPLYSNINGER – SKAL LÆSES INDEN BRUG  
BELANGRIJKE INFORMATIE – VOOR GEBRUIK DOORLEZEN  
OLULINE TEAVE – PALUN LUGEGE ENNE KASUTAMIST  
TÄRKEITÄ TIETOJA – LUETTAVA ENNEN KÄYTTÖÄ  
ΣΗΜΑΝΤΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ – ΠΑΡΑΚΑΛΟΥΜΕ ΔΙΑΒΑΣΤΕ ΠΡΙΝ ΑΠΟ ΤΗ ΧΡΗΣΗ  
FONTOS TUDNIVALÓK – HASZNÁLAT ELŐTT OLVASSA EL  
SVARĪGA INFORMĀCIJA – LŪDZU, IZLASIET PIRMS LIETOŠANAS  
SVARBI INFORMACIJA – PERSKAITYKITE PRIEŠ NAUDODAMI  
VIKTIG INFORMASJON – LES FØR BRUK  
WAŽNE INFORMACJE – NALEŻY PRZECZYTAĆ PRZED UŻYCIEM  
INFORMAÇÕES IMPORTANTES – LER ANTES DE UTILIZAR  
INFORMAȚII IMPORTANTE – VĂ RUGĂM SĂ LE CITIȚI ÎNAINTEA UTILIZĂRII  
DŮLEŽITÉ INFORMÁCIE – PREČÍTAJTE SI ICH PRED POUŽITÍM  
ПОМЕМБНЕ ИНФОРМАЦИЈЕ – ПРЕБЕРИТЕ ПРЕД УПОРАБО  
VIKTIG INFORMATION – LÄS FÖRE ANVÄNDNING  
ÖNEMLI BILGI – LÜTFEN KULLANMADAN ÖNCE OKUYUN

Contents (English - EN).....	1
Contenu (Français - FR).....	4
Contenuto (Italiano - IT) .....	8
Inhalt (Deutsch - DE) .....	12
Índice (Español - ES).....	16
内容 (日本語 - JA).....	20
Съдържание (Български - BG) .....	23
Sadržaj (Hrvatski - HR).....	27
Obsah (Česky - CS) .....	31
Indhold (Dansk - DA).....	35
Inhoud (Nederlands - NL).....	39
Sisukord (Eesti - ET).....	43
Sisältö (Suomi - FI) .....	47
Περιεχόμενα (Ελληνικά - EL).....	50
Tartalom (Magyar - HU) .....	54
Saturs (Latviski - LV) .....	58
Turinys (Lietuvių - LT).....	62
Innhold (Norsk - NO) .....	66
Spis treści (Polski - PL) .....	69
Índice (Português - PT).....	73
Cuprins (Română - RO).....	77
Obsah (Slovenčina - SK).....	81
Vsebina (Slovenščina - SL).....	85
Innehåll (Svenska - SV) .....	89
İçindekiler (Türkçe - TR).....	92

### **Instructions for Use/Intended Use**

The Tecomet hybrid cases and trays are intended for the protection, organization, and delivery to the surgical field of surgical instruments and/or other medical devices. The cases and trays are not intended to maintain sterility by themselves. They are designed to facilitate the sterilization process when used with an FDA approved sterilization wrap. Wrapping material is designed to allow air removal, steam penetration/evacuation (drying) and to maintain the sterility of the internal components.

### **Recommendations for Care, Cleaning, and Sterilization for the Acetabular Reamer Case**

Both physical and chemical (detergent) processes may be necessary to clean soiled items. Chemical (detergent) cleaners alone cannot remove all soil and debris; therefore, a careful manual cleaning of each item with a soft sponge or cloth is essential for maximum decontamination. For difficult areas, a clean soft bristled brush is recommended. Once the items have been cleaned, they should be thoroughly rinsed with clean water to remove any detergent or chemical residue before sterilization.

**Tecomet** recommends the use of a mild enzymatic detergent with a near neutral pH. Do not use solvents, abrasive cleaners metal brushes, or abrasive pads. Case and trays may be placed in mechanical cleaning equipment.

### **Materials and Restricted Substance**

For indication that the device contains a restricted substance or material of animal origin see product label.

**Always** inspect instruments and case for cleanliness and damage before use. Make sure all latches and handles are secure and in working order. Do not overload cases and trays. Balance contents uniformly within the container and arrange to allow steam to come into contact with all objects in the container. Tecomet recommends that the cases and trays be wrapped according to the sterilization wrap manufacturer's instructions prior to sterilization to maintain sterility of internal components/items and for proper aseptic presentation to the surgical field. Testing was performed with the Acetabular Reamer Case and the instruments included in the kit. Always follow instrument manufacturer instructions if their

sterilization or drying recommendations exceed these guidelines. For steam autoclaves, Tecomet has verified product performance in the following cycles.

### Pre-vacuum Sterilizer

Wrapped cases, trays, and instruments should be sterilized according to the table below.

Cycle Type	Temperature	Exposure Time	Dry Time	Cool-Down Time*
<b>United States Recommended Parameters</b>				
Pre-vacuum/ Vacuum Pulse	132°C / 270°F	4 minutes	40 minutes	30 minutes
Cycle Type	Temperature	Exposure Time	Dry Time	Cool-Down Time*
<b>European Recommended Parameters</b>				
Pre-vacuum/ Vacuum Pulse	134°C / 273°F	3 minutes	40 minutes	30 minutes

\*Allow the wrapped device to cool outside of the autoclave for 30 minutes prior to use.

### Filtered Sterilization Container Systems

Tecomet delivery systems have not been validated for use in filtered sterilization containers and Tecomet does not recommend use of such a system.

If filtered sterilization container systems are used, the user is responsible for following the recommendations of the manufacturer for the proper placement and use of cases and tray inside the container.

### Contraindications

The hybrid cases and trays have **NOT** been validated for use with flexible endoscopes, or devices with lumens or working channels longer than 4 inches (3mm ID). Always refer to instrument manufacturer instructions. The hybrid cases and trays have **NOT** been validated for ETO sterilization of devices. Always refer to instrument manufacturer instructions. Stacking the hybrid cases and overloading the cases will adversely affect sterilization and drying effectiveness. **DO NOT STACK** cases or trays in an autoclave chamber.

## **Notice**

1. **DO NOT** load cases into sterilizer on sides or upside down with lid side on shelf or cart. Load cases on cart or shelf so that the lid is always facing upwards. This will allow for proper drying. The hybrid cases and trays are designed to drain in this position.
2. **DO NOT STACK** cases or trays in an autoclave chamber.
3. After the autoclave door is opened, all cases must be allowed to cool thoroughly. Place cases on rack or shelf with linen cover until cooling is complete. The potential for condensation may increase if the case is not allowed to cool properly.
4. If condensation is observed, check to ensure steps 1, 2, and 3 were followed. In addition, verify that the steam, which is being used for sterilization processing, has a quality of more than 97%<sup>1</sup>. Also confirm that the sterilizers have been inspected for routine maintenance in accordance with manufacturer recommendations.
5. Silicone mat accessories should be placed in alignment with the perforations in the case or tray so that the sterilant pathways are not obstructed.

**The following table provides guidance on maximum loads. Do not overload delivery systems (cases and trays), and always follow AAMI, OSHA and hospital standards for maximum loading. Follow instrument manufacturer's instructions for use.**

<b>Tray Size</b>	<b>Maximum Total Weight</b>
Small (~ 4 x 8 x 1 in.)	1 lb. (0.45 kg)
Tray (~ 8 x 17 x 1 in.)	2 lb. (0.91 kg)
½ Size Case (~ 7 x 7 x 4 in.)	8 lb. (3.6 kg)
¾ Size Case (~ 7 x 12 x 4 in.)	13 lb. (5.9 kg)
Full Size Case (~ 9 x 19 x 4 in.)	22 lb. (10.0 kg)

## **REFERENCES**

<sup>1</sup>ANSI/AAIM ST79: Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities

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Symmetry Medical Manufacturing, Inc.
Symmetry Medical Polyvac S.A.S